



Dehider Sterilizer

PRODUCT DESCRIPTION

Industrial automatic sterilizer for dehiders used in meat harvest and fabrication. Commercial sterilization is achieved by spraying Dehider with 195 deg F hot water. Smart controls minimize water usage and ensures consistent results with all operators.

Includes slot to sterilize a knife.

Suitable for Jarvis or similar Dehiders.

BENEFITS

Dramatically reduces water usage - as little as 4.75-oz of water is used per sterilization cycle. Does not require a potable water rinse after sterilizing. When the Econoliser is not sterilizing, water is not being used.

Environmentally friendly - whether you use city or well water, reducing water usage minimizes waste and helps preserve the environment.

Reduces operating costs - includes cost reductions related to acquiring, pumping, heating, treating, and disposing water.

Performance - proven commercial sterilization. Dehider is sterilized with 195 deg F water. More effective in micro reduction than ozonated water or chemical sterilization.

Automatic - the sterilization cycle is activated by inserting Dehider. Achieves the same commercial sterilization results with all operators.



Safety - the hot water spray is contained within the sterilizer, removes the risk of scalding operators.

Reduces condensation - the hot water spray is enclosed and momentary, steam is not continuously released into surrounding environment.

Dehiders are left cleaner - the time in sterilization cycle and the pressurized hot water spray results in a cleaner power tool.



Dehider Sterilizer



SPECIFICATIONS

Temperature set point	194 deg F
Sterilization Cycle Time	4.5 sec (adj)
Water consumption - 4 nozzle unit	4.75 oz @ 4.5 sec cycle time

UTILITY REQUIREMENTS

Electrical connection	240/277 VAC, single phase, 60 Hz
Max power consumption	2.5 kVA or 2.5 kW
Max current	11 Amps
Water Connection	1/2-inch NPT
Water pressure (minimum)	45 psi
Water flow rate - 4 nozzle unit	1.0 gal/min